



Central pressure requires levelled samples

### CENTRAL PRESSURE

**SAMPLE HOLDER M (TEAR SHAPE)**  
for polishing head ECO



| made of stainless steel

Samples 5x Ø 20-40 mm  
Dimensions Ø 130 mm

Order No.: Z5445036

**SAMPLE HOLDER L (TEAR SHAPE)**

for polishing head ECO+,  
Qpol 250 BOT, Qpol XL,  
Saphir 550/560



| incl. adaption

Samples 6x Ø 20-40 mm  
Dimensions Ø 160 mm

Order No.: Z5446038

### CENTRAL PRESSURE

**SAMPLE HOLDER M (3-POINT-FIXATION)**

for polishing head ECO



| especially suitable for hot mounted samples  
| made of stainless steel

Size	Order No.:	Dimensions
8x Ø 25 mm	Z5445030	Ø 130 mm
8x Ø 30 mm	Z5445031	Ø 130 mm
5x Ø 40 mm	Z5445034	Ø 130 mm
8 x Ø 1"	Z5445030	Ø 130 mm
8 x Ø 1 1/4"	Z5445032	Ø 130 mm
6 x Ø 1 1/2"	Z5445033	Ø 130 mm
3x 55 x 30 mm	Z5445035	Ø 130 mm

**SAMPLE HOLDER L (3-POINT-FIXATION)**

for polishing head ECO+,  
Qpol 250 BOT, Qpol XL, Saphir 550/560



| incl. adaption

Size	Order No.:	Dimensions
12 x Ø 25 mm	Z5446030	Ø 160 mm
8 x Ø 30 mm	Z5446031	Ø 160 mm
6 x Ø 40 mm	Z5446034	Ø 160 mm
5 x Ø 50 mm	Z5446035	Ø 160 mm
3 x Ø 60 mm	Z5446042	Ø 185 mm
5 x Ø 70 mm	Z5446044	Ø 204 mm
12 x Ø 1"	Z5446030	Ø 160 mm
8 x Ø 1 1/4"	Z5446032	Ø 160 mm
6 x Ø 1 1/2"	Z5446033	Ø 160 mm
5 x Ø 2"	Z5446035	Ø 160 mm
3x 55 x 30 mm	Z5446036	Ø 160 mm
3x 60 x 40 mm	Z5446037	Ø 160 mm
3x 70 x 40 mm	Z5446043	Ø 185 mm
3x 100 x 50 mm	Z5446045	Ø 204 mm

**VARIO-SAMPLE HOLDER**



- | fixing, grinding and polishing of samples with big dimensions or complicated geometries
- | made of anodized aluminium
- | incl. adaption


**suitable sample holder on request**



Examples for application



# SPECIAL SAMPLE HOLDERS



The problem of clamping samples with special dimensions can be solved with special sample holders.

