

## Qprep GALAXY polishing cloths

QPREP GALAXY polishing cloths offer the right cloth material for all materials and analysis purposes. Low resilience is provided by our GALAXY carrier disc.



### PRODUCT ADVANTAGES

- | Faster application because there is no need to glue on the polishing cloth
- | Low resilience due to metallic carrier plate
- | Easier handling due to grip tabs on the carrier plate

### RECOMMENDED APPLICATIONS

- | Process-reliable and reproducible polishing results
- | Use with GALAXY magnetic foil



Item No.	Unit	Description
----------	------	-------------

### GALAXY POLISHING CLOTH WITH METAL BACK

#### Description Ø

#### GALAXY – ALPHA

- Suitable for pre-polishing
- Very hard, perforated chemical fiber cloth
- For ceramics, hard metal, steel, grey cast iron, aluminum
- High lifetime and stock removal
- For use with diamond grain size: 15/9 µm
- Fixation system: Magnetic foil

95001407	5 Pcs.	Alpha	200 mm
95001408	5 Pcs.	Alpha	250 mm
95001409	5 Pcs.	Alpha	300 mm
95013994	5 Pcs.	Alpha	350 mm

#### GALAXY – BETA

- Suitable for pre-polishing
- Hard synthetic cloth
- Materials with high hardness, steel, grey cast iron, hard metal, ceramics
- For use with diamond grain size: 15/9/6 µm
- Fixation system: Magnetic foil

95001410	5 Pcs.	Beta	200 mm
95001411	5 Pcs.	Beta	250 mm
95001412	5 Pcs.	Beta	300 mm
95013995	5 Pcs.	Beta	350 mm

#### GALAXY – DELTA

- Suitable for pre- & intermediate polishing
- Medium hard silk cloth
- C-steels, precious metal, plastics [CFK, GFK], coatings, aluminum
- For use with diamond grain size: 9/6/3 µm
- Fixation system: Magnetic foil

95001413	5 Pcs.	Delta	200 mm
95001414	5 Pcs.	Delta	250 mm
95001415	5 Pcs.	Delta	300 mm
95013996	5 Pcs.	Delta	350 mm