

| Item No. | Unit | Description | |
|--|--------|-------------|--------|
| GALAXY POLISHING CLOTH WITH METAL BACK | | | |
| Description Ø | | | |
| GALAXY – PHI <ul style="list-style-type: none"> • Suitable for pre-, intermediate-, and end polishing • Medium hard chemical fiber cloth • For iron, non – ferrous metals, minerals, ceramics, composite materials • For use with diamond grain size: 9/6/3/1 µm • Fixation system: Magnetic foil | | | |
| 95012739 | 5 Pcs. | Phi | 200 mm |
| 95012740 | 5 Pcs. | Phi | 250 mm |
| 95012741 | 5 Pcs. | Phi | 300 mm |
| 95014002 | 5 Pcs. | Phi | 350 mm |
| GALAXY – SIGMA <ul style="list-style-type: none"> • Suitable for intermediate and end polishing • Soft wool cloth • For all materials • For use with diamond grain size: 6/3/1 µm • Fixation system: Magnetic foil | | | |
| 95001416 | 5 Pcs. | Sigma | 200 mm |
| 95001417 | 5 Pcs. | Sigma | 250 mm |
| 95001418 | 5 Pcs. | Sigma | 300 mm |
| 95013997 | 5 Pcs. | Sigma | 350 mm |
| GALAXY – ETA <ul style="list-style-type: none"> • Suitable for end polishing • Short flocked, medium hard synthetic cloth • For all materials of medium and high hardness • For use with diamond grain size: 3/1 µm • Fixation system: Magnetic foil | | | |
| 95016832 | 5 Pcs. | Eta | 200 mm |
| 95016833 | 5 Pcs. | Eta | 250 mm |
| 95016834 | 5 Pcs. | Eta | 300 mm |
| 95016835 | 5 Pcs. | Eta | 350 mm |
| GALAXY – IOTA <ul style="list-style-type: none"> • Suitable for end polishing as well as fine polishing with fine polishing suspension • Short flocked, soft synthetic cloth • For all materials, especially hard materials • For use with diamond grain size: 3/1 µm • Suitable for oxide suspensions • Fixation system: Magnetic foil | | | |
| 95001419 | 5 Pcs. | Iota | 200 mm |
| 95001420 | 5 Pcs. | Iota | 250 mm |
| 95001421 | 5 Pcs. | Iota | 300 mm |
| 95013998 | 5 Pcs. | Iota | 350 mm |
| GALAXY – ZETA <ul style="list-style-type: none"> • Suitable for end polishing as well as fine polishing with fine polishing suspension • Short flocked, soft synthetic cloth • For all materials • For use with diamond grain size: 3/1 µm • Suitable for oxide suspensions as well as alumina • Fixation system: Magnetic foil | | | |
| 95005776 | 5 Pcs. | Zeta | 200 mm |
| 95005777 | 5 Pcs. | Zeta | 250 mm |
| 95005778 | 5 Pcs. | Zeta | 300 mm |
| 95013999 | 5 Pcs. | Zeta | 350 mm |

