

Aprep Polishing cloths, self-adhesive

QPREP polishing cloths offer the right cloth material for all materials and analysis purposes.



PRODUCT ADVANTAGES

- I Ideal for polishing in combination with the media carrier for the fully automatic Qpol 250 BOT as well as with the X-Tap system
- I Optimized adhesive backing ensures uniform fixing over the entire surface on the fixation system and realized no slipping or detachment of the polishing cloth

RECOMMENDED APPLICATIONS

- For use with X-Tap system
- I For use with media carrier plate (Qpol 250 BOT)

	ltem No.	Unit	Description POLISHING CLOTHS, SELF-ADHESIVE Description Ø ALPHA – Polishing cloth • Suitable for pre-polishing • Very hard, perforated chemical fiber cloth • For ceramics, hard metal, steel, grey cast iron, aluminum • High lifetime and stock removal • For use with diamond grain size: 15/9 μm • Fixation system: X-Tap		
	92002555	5 Pcs.	Alpha	200 mm	
	92002564	5 Pcs.	Alpha	250 mm	
	92002573	5 Pcs.	Alpha	300 mm	
	95016029	5 Pcs.	Alpha	350 mm	
			ceramics • For use with diamond grain size: 15/9/6 µm • Fixation system: X-Tap		
	95006571	5 Pcs.	Beta	200 mm	
	95006572	5 Pcs.	Beta	250 mm	
	95006573	5 Pcs.	Beta	300 mm	
	95016030	5 Pcs.	Beta	350 mm	
3 3 9			DELTA – Polishing cloth • Suitable for pre- & intermediate polishing • Medium hard silk cloth • C-steels, precious metal, plastics [CFK, GFK], coatings, aluminum • For use with diamond grain size: 9/6/3 µm • Fixation system: X-Tap		
	92008801	5 Pcs.	Delta	200 mm	
	92008802	5 Pcs.	Delta	250 mm	
	92008803	5 Pcs.	Delta	300 mm	