

Item No.	Unit	Description	
POLISHING CLOTHS, SELF-ADHESIVE			
Description Ø			
GAMMA – Polishing cloth <ul style="list-style-type: none"> • Suitable for pre- & intermediate polishing • Medium hard silk cloth • C-steels, precious metal, coatings, plastics [CFK, GFK], • For use with diamond grain size: 9/6/3 µm • Fixation system: X-Tap 			
92002558	5 Pcs.	Gamma	200 mm
92002567	5 Pcs.	Gamma	250 mm
92002576	5 Pcs.	Gamma	300 mm
95016032	5 Pcs.	Gamma	350 mm
PHI – Polishing cloth <ul style="list-style-type: none"> • Suitable for pre-, intermediate-, and end polishing • Medium hard chemical fiber cloth • For iron, non – ferrous metals, minerals, ceramics, composite materials • For use with diamond grain size: 9/6/3/1 µm • Fixation system: X-Tap 			
92005678	5 Pcs.	Phi	200 mm
92005679	5 Pcs.	Phi	250 mm
92005680	5 Pcs.	Phi	300 mm
95016031	5 Pcs.	Phi	350 mm
SIGMA – Polishing cloth <ul style="list-style-type: none"> • Suitable for pre- & intermediate polishing • Soft wool cloth • For all materials • For use with diamond grain size: 6/3/1 µm • Fixation system: X-Tap 			
92008810	5 Pcs.	Sigma	200 mm
92008811	5 Pcs.	Sigma	250 mm
92008812	5 Pcs.	Sigma	300 mm
95016033	5 Pcs.	Sigma	350 mm
ETA – Polishing cloth <ul style="list-style-type: none"> • Suitable for end polishing • Short flocked, medium hard synthetic cloth • For all materials of medium and high hardness • For use with diamond grain size: 3/1 µm • Fixation system: X-Tap 			
95016828	5 Pcs.	Eta	200 mm
95016829	5 Pcs.	Eta	250 mm
95016830	5 Pcs.	Eta	300 mm
95016831	5 Pcs.	Eta	350 mm
IOTA – Polishing cloth <ul style="list-style-type: none"> • Suitable for end polishing as well as fine polishing with fine polishing suspension • Short flocked, soft synthetic cloth • for all materials, especially hard materials • For use with diamond grain size: 3/1 µm • Suitable for oxide suspensions • Fixation system: X-Tap 			
95002393	5 Pcs.	Iota	200 mm
95002394	5 Pcs.	Iota	250 mm
95002395	5 Pcs.	Iota	300 mm
95016036	5 Pcs.	Iota	350 mm

