



Item No.	Unit	Description
<div style="display: flex; align-items: center;"> <div style="margin-right: 10px;"> <div style="border: 1px solid black; padding: 2px; display: inline-block;">5 Units / -5%</div> <div style="border: 1px solid black; padding: 2px; display: inline-block; margin-left: 5px;">10 Units / -10%</div> </div> <div style="background-color: #800000; color: white; padding: 5px; border-radius: 5px;">Hot mounting material</div> </div>		
<div style="display: flex; align-items: center;"> <div style="margin-right: 10px;"> <div style="background-color: red; color: white; border-radius: 50%; padding: 5px; display: inline-block; font-weight: bold;">NEW</div> </div> <div> <p>EPO-Max black, epoxy resin for gap-free embedding » very high hardness, easy cleaning of mould and ram due to low adhesion » with mineral filler</p> </div> </div>		
95013811	1 kg	EPO-Max black
95013812	5 kg	EPO-Max black
95013813	10 kg	EPO-Max black
<p>Epo black, epoxy resin, suppressed gap forming » very high hardness » contains a high proportion of glass and mineral filler</p>		
95011990	1 kg	Epo black
95011991	5 kg	Epo black
95011992	10 kg	Epo black
<p>Duroplast blue, diallylphthalate, suppressed gap forming » general purpose and especially for grinding stones » high hardness » contains a high proportion of glass fiber</p>		
95011999	1 kg	Duroplast blue
95012000	5 kg	Duroplast blue
95012001	10 kg	Duroplast blue
<p>Duroplast black, phenolic resin for SEM examinations and electrolytic polishing » medium hardness » conductive, contains graphite particles</p>		
95011993	1 kg	Duroplast black
95011994	5 kg	Duroplast black
95011995	10 kg	Duroplast black
<p>Thermoplast, acrylic resin for transparent mounts » ideal for targeted preparations » low - medium hardness</p>		
95011996	1 kg	Thermoplast
95011997	5 kg	Thermoplast
95011998	10 kg	Thermoplast
<p>Bakelite black, phenolic resin for general purpose » medium hardness » contains wood powder » very convenient material for filling</p>		
95011981	1 kg	Bakelit black
95011982	5 kg	Bakeliet black
95011983	10 kg	Bakelite black
<p>Bakelite red, phenolic resin for general purpose » medium hardness » contains wood powder</p>		
95011984	1 kg	Bakelite red
95011985	5 kg	Bakelite red
95011986	10 kg	Bakelite red
<p>Bakelite green, phenolic resin for general purpose » medium hardness » contains wood powder</p>		
95011987	1 kg	Bakelite green
95011988	5 kg	Bakelite green
95011989	10 kg	Bakelite green

