

Data sheet

Laboratory Shakers LS 500



The photo above is for reference only, may show additional options not included in standard equipment. The real appearance, particularly color and structure of the material may differ from the ones presented in the photo.

Standard equipment:

- microprocessor control of rotation and time
- LCD digital display
- virable speed control: 30 ... 300rpm
- shaking mode: from 1min to 99h or continuous operation

TECHNICAL DATA	
movement	orbital
controller	microprocessor
display	LCD
speed range [rpm]	30300
acurracy [rpm]	10
amplitude [mm]	5 or 10
max load capacity[kg]	10
shaking mode	1min 99h or continuous operation
dimensions without/with shaking table [mm]	
width	550
heigth	120/220
depth	440
fits to cooled incubator	ILW 240
nominal power [W]	60
weight with shaking table [kg]	22
ambient temperature [°C]	+10+40
humidity [%]	up to 70
voltage	230 V 50 Hz
warranty	24 months
manufacturrer	POL-EKO-APARATURA



ACCESSORIES	
universal platform	
platform for Petri plates shaking	
platform for separatory funnels	
platform for fixing flasks handles	
anti-skid mat	