

Data sheet

Cooled Incubator (ST) ST 2 Smart PRO



The photo above is for reference only, may show additional options not included in standard equipment.

The real appearance, particularly color and structure of the material may differ from the ones presented in the photo


Advantages of the SMART PRO controller:

- large (7"), clear, full colour touch screen
- LAN, USB ports and WiFi for communication and data transfer
- multi-segment time and temperature programs
- overview of data in tabular and graphic form
- visual and sound alarm
- Admin function for management
- password protected log-in
- internal memory for programs and data storage
- event registry with user notifications
- LabDesk software and user manual for direct download

TECHNICAL DATA		
air convection	forced	
chamber capacity [I]	150	
working capacity [I]	122	
controller	microprocessor with external LCD graphic display	
display	7" full colour touch screen	
TEMPERATURE		
temperature range [°C]	+3+70	
temperature resolution every [°C]	0,1	
over temperature protection	class 3.3 to DIN 12880	



CHAMBER	
door type	solid / glass or double (option) /1/
interior	
P Smart PRO	acid-proof stainless steel to DIN 1.4301
PS Smart PRO	acid-proof stainless steel to DIN 1.4301
housing	
P Smart PRO	powder coated sheet
PS Smart PRO	stainless steel polished
overall dims [mm] /2/	
width A	620
height B	900
depth C	650
internal dims [mm] /3/	
width D	480
width D'	520
height E	660
depth F	420
depth F'	480
depth G	320
height H	440
height I	-



shelves (standard / max)	3 // 4	
max shelf workload [kg] /4/	10	_
- reinforced shelf version (PW) /5/	-	
max unit workload [kg]	30	
- reinforced unit version (W) /6/	-	
weight [kg]	54	

230V 50Hz
250
24 months
POL-EKO-APARATURA
_

all the above technical data refer to standard units (without optional accessories) * - 230V 60Hz, 115V 60Hz also available

- additional internal glass door
 depth doesn't include 50 mm of power cable
 dims of units with double doorcan be smaller
- 4 on uniformly loaded surface 5 reinforced shelf
- 6 reinforced version



OPTIONS AND ACCESSORIES



Order number: */C Internal glass door



Order number: */A External glass door



Order number: */P INOX Stainless steel wire shelf



Order number: */PP Perforated shelf



Order number: */PW Reinforced shelf



Order number: KUW GN*/* Stainless steel cuvettes



Order number: ST/CHL/SWP ALU Aluminum drawer with powder coated slides



number: ST/CHL/SWP INOX Stainless steel drawer with powder coated slides



Order number: ST/CHL/SWPN INOX Stainless steel drawer with stainless steel slides



Order number: ORG-FARM Pharma organizer



Order number: QLK* Wheels



Order number: */S or */S INOX Table with wheels



Order number: */ST or */S INOX Table with wheels



Order number: *PLUS Automatic defrosting function



Order number: PT 100 Additional temperature sensor



Order number: GNZ Internal socket



Order number: PHR Humidity measurement



Order number: OCZ/N Non-standard access port for external sensor



Order numberLabDesk LabDesk - software



Order number: BRT/*/L or IQ/OQ/PQ Calibration and IQ, OQ, PQ qualification